

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## FLAER: Alexa fluor<sup>®</sup> 488 proaerolysin

**FL2 – 50 µg lyophilized FLAER; reconstitute in 1.0 ml of PBS**

### **Description:**

The reagent is now supplied as a powder, which after reconstitution, produces a  $10^{-6}$  M FLAER stock solution. This solution also contains 5 mg/mL albumin, added in order to maximize stability. Vials of FLAER should be kept in a freezer protected from light before use. Sodium azide may be added as a preservative. **The stock solution may be kept in the refrigerator protected from light for at least 1 year.**

FLAER is produced by labeling an inactive variant of the protein proaerolysin with Alexafluor<sup>®</sup> 488 succinimidyl ester, provided under agreement with Life Technologies Corporation. The reagent binds selectively to cell surface GPI-anchored proteins. It is intended for research purposes only and is not for resale.

$5 \times 10^{-8}$  M has proven to be a useful final concentration of FLAER for flow cytometry (this represents a 1 to 20 dilution of the stock directly into the test sample), however you may wish to try other concentrations for your own application.

**There are a number of published methods describing the detection of PNH using FLAER. One example is Sutherland et al. in Cytometry Part B (Clinical Cytometry) 72B: 167-177 (2007). We can provide a detailed protocol upon request.**

Manufactured by:

**PINEWOOD SCIENTIFIC SERVICES INC**  
210 Denison Rd., Victoria, BC, V8S 4K3, Canada  
Tel: 250-598-6822, Fax: 250-598-6877  
protoxbiotech.com

Updated: January 19, 2011

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)